


GSV – GERICKE SAMPLING VALVE

Representative and automated sampling high potent OSD production

In the pharmaceutical industry, ensuring the quality of oral solid dosage (OSD) forms is paramount. Gericke innovative sampling valve plays a crucial role in continuous manufacturing processes, providing precise and representative quality assurance. It works in a fully automated fashion and is the only sampler on the market which takes the complete cross-sectional probe to provide a representative sample. It is designed to meet the rigorous demands of pharmaceutical production. It allows for accurate sampling of small material quantities during the manufacturing process, ensuring that each batch meets the highest quality standards.

GERICKE AUTOMATED SAMPLER – GSV 65 CA



ADVANTAGES

- **Precision:** Accurate and representative sampling ensures output meets quality standards before further processing.
- **Integration:** Seamlessly integrates with continuous manufacturing systems, allowing automated sampling without production interruptions.
- **Compliance:** Fully cGMP compliant designed allowing for fast cleaning and product change over.
- **Efficiency:** Adjustable sampling size to meet unit dose, minimizing material waste.
- **Versatility:** Suitable for various oral solid dosage forms, providing a simple, safe, and effective solution.
- **Safety:** Suitable for high containment applications with Wash-in-Place option for safe dismantling.